



Your Well-being Companion

www.companio.in

ISO 13485 Certified Company

FOOT PULSE EXO



FOOT PULSE EXO

Foot Pulse EXO is a 30 minutes auto-programmed electrical impulse therapy. It is based on the combination of 12 different Electrical Muscle Stimulation programs with up to 99 level of intensity.

How Does It Work

In this therapy, the mild electrical pulses are sent through the soles of the feet. This causes the muscles in the lower leg to expand and contract, pushing the blood from the leg to the heart. The constant and regular stimulation boosts the blood circulation and helps open blockages in the arteries and veins. The increased blood circulation ensures that adequate oxygen and nutrients reach the muscles and help to heal pain, strengthen the muscles and make them flexible.



Foot Pulse ACU



Benefits Of Foot Pulse ACU



- Relief in Swollen feet and Ankles



- Helps in Muscle Strengthening



- Muscle Rehabilitation



- Relief in Muscle Cramps in Legs



- Relief in Calf Cramps

Know Your Product



- ① Sole / Foot Pads
- ② Time
- ③ Foot Intensity
- ④ Body Patches Intensity
- ⑤ Time Buttons
-/+ 1-60 minutes
(default 30 minutes)
- ⑥ Sole Intensity Level Buttons
-/+ 1-99 (Electrode pads)
- ⑦ Body Intensity Level Buttons
-/+ 1-99 (Electrode pads)
- ⑧ Power Button
Press to turn device on/off
- ⑨ Power and Accessory Sockets
Power and Electrode Pads sockets (simply plug in the leads, connect the pads and select body intensity levels)
- ⑩ AC/DC Adaptor
- ⑪ Remote Control
- ⑫ Electrode Pads
To target specific areas of pain and swelling

Relief From Various Body Issues

companio[®]
Your Well-being Companion



Electrode Body Patches

Use It With
Foot Pulse Body Function



Back Pain



Calf Pain




Thigh Pain



Neck Pain




Shoulder Pain



Who Cannot Use This Product




Pregnant Women



If You Have A Pacemaker



Epileptic People



People With Deep Vein Thrombosis



Those With A Heart Condition

Accessories



Electrode Gloves



Electrode Gel




Multi
Electrode Patches



Bag



How To Use ?



STEP#1
Plug the power cable into the nearest electricity socket and the power input of the device.




STEP#2
Sit down on a chair or a sofa and place your both bare feet on the foot shaped pad of the device.



STEP#3
To turn the device ON, press the power button on the device or the remote control. Adjust the intensity through -/+ button.



STEP#4
Connect the electrode pads with the connecting wire. Now plug this wire into the accessories socket on the device.



STEP#5
Peel off the protective membrane from the body pads and carefully place them on the clean, dry, uninjured skin.

Certification



Manufacturing Unit





Your Well-being Companion

Manufacturer

CAMEX WELLNESS LIMITED

16 &17, Akshar Industrial Park, Nr. Vasna - Chacharwadi Village,

Changodar - Bavla Highway, Changodar, Sanand - 382110

+91 9898852989 | welwisher@camexwellness.com | www.companio.in