



Your Well-being Companion

www.companio.in

ISO 13485 Certified Company

BACK PULSE




BACK PULSE

Back Pulse is a 18 minutes auto-programmed electrical impulse therapy. It is based on the combination of 7 different Transcutaneous Electrical Nerve Stimulation programs with up to 20 level of intensity.

How Does It Work

The TENS technology delivers mild electrical impulses to nerve fibers and thereby blocks the pain messages from reaching the brain. Also, they produce natural pain-killing hormones Endorphins that effectively alleviates the back pain and make you feel energized and active.




Anytime, Anywhere Access


Simple To
Use



Easy To
Wear



No Side
Effects



Benefits Of Back Pulse



- Reduce Back Pain




- Build Back Strength




- Improve Back Movement

Know Your Product



Relief From Various Back Issues


companio[®]
Your Well-being Companion



Who Cannot Use This Product




Pregnant Women




If You Have A Pacemaker




Epileptic People




How To Use ?



STEP#1
Lay the belt.




STEP#2
Remove the replaceable electrode gel from the foil.




STEP#3
Stick the replaceable electrode gel.



STEP#4
Wear the belt on the lower back.



STEP#5
Press switch on button to start.



STEP#6
Press '+' to increase and '-' to decrease the intensity.

Certification

