



Your Well-being Companion

www.companio.in

ISO 13485 Certified Company

KNEE PULSE



KNEE PULSE

Knee Pulse is a 15 minutes auto-programmed electrical impulse therapy. It is based on the combination of 7 different Transcutaneous Electrical Nerve Stimulation programs with up to 20 level of intensity.


How Does It Work



The device releases therapeutic impulses that work to suppress the transmission of pain in the nerve fibers of the knee muscles and also release the pain-relieving hormones - Endorphins to alleviate pain. ACU technology is a tested, effective and medication-free method to treat pain without any harmful side-effects.

Anytime, Anywhere Access

Simple To
Use



Easy To
Wear



No Side
Effects



Benefits Of Knee Pulse



- Reduces knee osteoarthritis Pain
- Reduces inflammation of Knees
- Helps increase muscle Strength
- Boosts blood circulation and Stimulates nerves around Knee
- Helps in slowing down or Preventing joint damage With pain relief




Know Your Product



- ① Operating Unit
- ② On / Off
- ③ Selection key - Increase Intensity (+) Plus sign
- ④ Selection key - Reduce Intensity (-) Minus sign
- ⑤ Pads
- ⑥ Electrode Gel and Replacement

Relief From Various Knee Issues

companio[®]
Your Well-being Companion



Who Cannot Use This Product




Pregnant Women




Epileptic People


People With Deep Vein Thrombosis




On The Region Of Skin Injuries, Swelling, Bruise And Cuts




How To Use ?



STEP#1
Lay the Knee Pulse belt.



STEP#2
Remove the replaceable electrode gel from the foil.



STEP#3
Stick the replaceable electrode gel.



STEP#4
Wear the belt on your knee and make sure the patella is in between.



STEP#5
Press switch on button to start.




STEP#6
Press '+' to increase and '-' to decrease the intensity.

Certification



Manufacturing Unit





Your Well-being Companion

Manufacturer

CAMEX WELLNESS LIMITED

16 &17, Akshar Industrial Park, Nr. Vasna - Chacharwadi Village,

Changodar - Bavla Highway, Changodar, Sanand - 382110

+91 9898852989 | welwisher@camexwellness.com | www.companio.in